

TECHNICAL INFORMATION SHEET

SOCLEAN ANTIVIRAL+ SANITISER

Product Description

A low toxicity dual purpose product which cleans and sanitises all hard surfaces such as floors, paintwork, working surfaces, stainless steel etc. in one operation.

Features and Benefits

- An excellent general-purpose sanitiser based on a natural organic biocide and plants extracts.
- Enables cleaning and sanitising to be carried out simultaneously.
- Effective in all water conditions.
- Can be used for cleaning and sanitising food utensils, glasses, crockery and all types of kitchen equipment.
- Suitable for use on all types of hard surfaces including plastics, laminates, glass and stainless steel.
- No residual organoleptic effect when used to clean stainless steel surfaces under normal cleaning conditions.
- Efficiently hinders lime scale build-up.
- Suitable for use in high risk infection areas e.g. public washrooms, healthcare environments.
- Tested, according to the most recent European Standards (EN-Norms) and proven to be an outstanding disinfectant with bactericidal and virucidal properties.

Use Instructions

General cleaning and disinfection to EN1276.

1. Use neat.
2. For heavily soiled areas, pre clean surfaces.
3. Apply via trigger spray, cloth or brush.
4. Leave for at least 30 seconds and wipe with a clean cloth.

Safe Handling and Storage

See visual information to right.

Full guidance on the handling and safe use of the product can be obtained from the Safety Data Sheet.

Store in original container.

Avoid extremes of temperature.

Health & Safety

HAZARD AS SUPPLIED



Not classified. Wear gloves for prolonged use.



HAZARD WHEN USED AT RECOMMENDED DILUTION.



Not classified. Wear gloves for prolonged use.



Technical Information

Appearance: Pale brown liquid

pH : <5.0

Relative density: 1.01

Odour: Characteristic

Solubility: Fully miscible in water

The above data is typical of normal production and should not be taken as a specification.

Materials Compatibility

Under recommended conditions of use suitable for use on most materials commonly encountered. Do not use on marble or other acid sensitive materials.

Other Information

Do not mix with conventional detergents as they may reduce product effectiveness.

Efficacy Data

Tested independently in accordance with:

- EN 1276 Bacterial action against *Pseudomonas aeruginosa* , *Escherichia coli*, *Staphylococcus aureus*, *Enterococcus hirae*, 30 seconds, dirty conditions.
- EN 13704 Sporicidal activity against *Clostridium difficile*, 5 minutes, clean conditions.
- EN14476 Virucidal activity against Norovirus, Hepatitis C, 5 minutes, clean conditions.
- EN14476 Virucidal activity against HIV, 15 minutes, clean conditions.
- EN13727 Bacterial action against the NDM1 bacterial strain, 1 minute, dirty conditions.
- TES-S-004: Triangle test for potential taint due to contact with materials.

A summary of the tests is available on request.